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TABLE 4-V. Character Table for C_{3v} Group

	$E(1)$	$C_3(2)$	$\sigma_v(3)$
a_1	1	1	1
a_2	1	1	-1
e	2	-1	0
Variables and Characters for Calculating Allowed Infrared Bands			
ϕ	0°	120°	0°
χ_m	3	0	1
U_R	7	4	5
θ	15	0	5

$$N_{a1} = \frac{1}{6}[1 \cdot 3 \cdot 1 + 2 \cdot 0 \cdot 1 + 3 \cdot 1 \cdot 1] = 1 \quad N_{a2} = \frac{1}{6}[1 \cdot 3 \cdot 1 + 2 \cdot 0 \cdot 1 + 3 \cdot 1 \cdot (-1)] = 0$$

$$N_e = \frac{1}{6}[1 \cdot 3 \cdot 2 + 2 \cdot 0 \cdot (-1) + 3 \cdot 1 \cdot 0] = 1$$

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TABLE 4-VI. Selection Rules for Binary Combination Frequencies of a Molecule with C_{3v} Symmetry

Combination	Number of times the vibration type appears in the character of the combination			Activity*	
	a_1	a_2	e	Raman	Infrared
$a_1 \times a_1$	1	0	0	+	+
$a_1 \times a_2$	0	1	0	-	-
$a_1 \times e$	0	0	1	+	+
$a_2 \times a_2$	1	0	0	+	+
$a_2 \times e$	0	0	1	+	+
$e \times e$	1	1	1	+	+

* + indicates active; - indicates inactive.

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$$N = \frac{1}{N_g} \sum N_e \theta X_i$$

Why FT-IR?

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
Fourier Transform Infrared (FT-IR) spectrometry was developed in order to overcome the limitations encountered with dispersive instruments. The main difficulty was the slow scanning process. A method for measuring all of the infrared frequencies **simultaneously**, rather than individually, was needed. A solution was developed which employed a very simple optical device called an **interferometer**. The interferometer produces a unique type of signal which has all of the infrared frequencies “encoded” into it. The signal can be measured very quickly, usually on the order of **one second** or so. Thus, the time element per sample is reduced to a matter of a few seconds rather than several minutes.

Most interferometers employ a **beamsplitter** which takes the incoming infrared beam and divides it into two optical beams. One beam reflects off of a flat mirror which is fixed in place. The other beam reflects off of a flat mirror which is on a mechanism which allows this mirror to move a very short distance (typically a few millimeters) away from the beamsplitter. The two beams reflect off of their respective mirrors and are recombined when they meet back at the beamsplitter. Because the path that one beam travels is a fixed length and the other is constantly changing as its mirror moves, the signal which exits the interferometer is the result of these two beams “interfering” with each other. The resulting signal is called an **interferogram** which has the unique property that every data point (a function of the moving mirror position) which makes up the signal has information about every infrared frequency which comes from the source.

This means that as the interferogram is measured, all frequencies are being measured **simultaneously**. Thus, the use of the interferometer results in extremely fast measurements.


Because the analyst requires a **frequency spectrum** (a plot of the intensity at each individual frequency) in order to make an identification, the measured interferogram signal can not be interpreted directly. A means of “decoding” the individual frequencies is required. This can be accomplished via a well-known mathematical technique called the **Fourier transformation**. This transformation is performed by the computer which then presents the user with the desired spectral information for analysis.

The Sample Analysis Process


The normal instrumental process is as follows:

- 1. The Source:** Infrared energy is emitted from a glowing black-body source. This beam passes through an aperture which controls the amount of energy presented to the sample (and, ultimately, to the detector).
- 2. The Interferometer:** The beam enters the interferometer where the “spectral encoding” takes place. The resulting interferogram signal then exits the interferometer.
- 3. The Sample:** The beam enters the sample compartment where it is transmitted through or reflected off of the surface of the sample, depending on the type of analysis being accomplished. This is where specific frequencies of energy, which are uniquely characteristic of the sample, are absorbed.
- 4. The Detector:** The beam finally passes to the detector for final measurement. The detectors used are specially designed to measure the special interferogram signal.
- 5. The Computer:** The measured signal is digitized and sent to the computer where the Fourier transformation takes place. The final infrared spectrum is then presented to the user for interpretation and any further manipulation.

Because there needs to be a relative scale for the absorption intensity, a **background** spectrum must also be measured. This is normally a measurement with no sample in the beam. This can be compared to the measurement with the sample in the beam to determine the “percent transmittance.” This technique results in a spectrum which has all of the instrumental characteristics removed. Thus, all spectral features which are present are strictly due to the sample. A single background measurement can be used for many sample measurements because this spectrum is characteristic of the **instrument** itself.

A Simple Spectrometer Layout

How does FTIR work? (www.ijvs.com)

The basic idea

The building blocks of an F-T are very different from the predecessor, the double beam (diffraction) infrared spectrometer. In an F-T you have the chain –

Source → Interferometer → Sample → Detector

Each building block varies from manufacturer to manufacturer but basically they all obey the same laws of physics. I can't be exhaustive, but I will try to cover most of the bits and pieces used in contemporary instruments.

Sources

The easy bit. All the manufacturers use a heated ceramic source. The composition of the ceramic and the method of heating vary but the idea is always the same, the production of a heated emitter operating at as high a temperature as is consistent with a very long life. Historically, people have achieved this objective in almost bizarre ways e.g. using an old fashioned gas mantle, or the Nernst Glower [a device utilising the fact that some ceramics conduct electricity when heated. The Nernst had to be pre-heated and then could be used as a conductor and hence as a hot source – they were tricky to use and unreliable]. These days, the manufacturers tend to go for either a conducting ceramic or a wire heater coated with ceramic. Now, as I explained recently when talking about hot samples in Raman spectroscopy, heated objects emit, the emission occurs at all wavelengths AND the emission at any wavelength increases with temperature. The plots of emission vs. frequency are shown in Figure 1.


Figure 1. Blackbody curves for a variety of temperatures on an absolute wavenumber abscissa scale, (a) in watts, (b) in counts per second. As the temperature is increased the overall intensity increases and the maximum of the curve shifts to higher temperature.

All infrared spectrometers are energy limited i.e. they are short of photons hence it is obvious that the hotter the source, the better. Equally obviously, it would be disastrous if the source operated at a variable temperature so manufacturers go to some lengths to prevent draughts passing the hot source. This is not as easy

as it sounds because air will inevitably convect past the source due to its temperature. Many manufacturers get round this problem by monitoring the source output, using the output signal to control the electrical power in a feedback loop.


Figure 2. The Source emits and is detected by D_s . Output regulates the chip C_D which controls the power fed to the source.

Ideally, then we have a stabilised hot source operating at a feverishly high temperature. Well not quite!

As all you folks who have kept an eye on the atmospheric pollution game will know, heating air to high temperatures can be problematic – generation of the notorious $(NO)_x$!

The equilibrium $N_2 + O_2 \rightleftharpoons NO_2 \rightleftharpoons N_2O_4$ lies overwhelmingly towards the reactants BUT moves towards NO_2 and hence N_2O_4 at elevated temperatures. Unfortunately, as it turns out, good modern sources can be hot enough to generate measurable amounts of NO_2 (and other nitrogen oxides hence $(NO)_x$). In the atmosphere, this would be only a trivial nuisance, one swamped by other sources of pollution but inside the sealed housing of an FTIR, the $(NO)_x$ can well cause unwanted absorption and even corrosion of the metal and plastic components. So, sources are hot but not too hot. Temperatures vary but a little above $1000^\circ C$ are fairly typical. Talk of $2000K$ is, I suspect, not realistic.

Interferometers

We all know that the interferometer is the heart of our instrument. The interferometer is the bit that analyses the infrared or near infrared and hence enables us to generate a spectrum. Most of us will have at one time or another seen the classical diagram of an interferometer and had it explained to us. Let's go through this again.

The classic Michelson Interferometer involves a beam splitter – a component which reflects about $\frac{1}{2}$ of the radiation that hits it and transmits the rest. One bundle of radiation follows one path and the remainder a different path. The two bundles are then recombined and pass out of the device. Now, one path is changed in length allowing radiation to interfere with that following the fixed path. In Figure 3, I show the classical diagram of the Michelson Interferometer.


Figure 3.

Parallel radiation entering the device hits beam splitter BS and is either reflected to mirror M_F or transmitted to hit mirror M_M . Radiation reflected by the mirrors then passes either through or is reflected off the beam-splitter and hence either goes back towards the source OR goes forward and out of the interferometer. Thus, only half the radiation entering the interferometer can get out in the OUT direction. Half is reflected back towards the input.

Let's assume we use monochromatic radiation of wavelength λ . If the paths involving M_F and M_M are identical in length, about $\frac{1}{2}$ of the radiation will pass through the interferometer, but if the paths differ by $\lambda/2$ when recombined on the beam-splitter, the two bundles of radiation will be out-of-phase and hence cancel each other out. Now, if M_M is moved, the efficiency with which the interferometer will pass radiation follows what is shown in Figure 4.


Figure 4.

If now, we use radiation of 2 wavelengths λ_1 and λ_2 the result is shown in Figure 5,


Figure 5.

that is, we will actually observe the sum of the two cosine functions. If white light is used, we will see the sum of all the cosine functions typical of ALL the wavelengths.

Now, note one feature – at zero path difference all the radiation of whatever wavelength passes through the interferometer. At any other path difference, only some will pass. Hence, the plot of the efficiency of transmissions vs. path difference for white light looks something like



Figure 6.

And is characterised by a "centre burst" at zero path difference and a very complex pattern of waves symmetrically dispersed about it. [Go now to your FTIR and call up the interferogram of the background and you will see a much better picture than the one I have drawn for you in Figure 6].

The interferogram, as Figure 6 is not recognisable as a spectrum. To convert it into something useful, it needs to be **Fourier Transformed**. At this stage, all textbooks launch into calculus, but I have always found that students went cross-eyed when I tried, so here is a simple-minded explanation of the Fourier Transform process.

Let's assume the interferogram is digitised and stored in a processor. Now let us also assume that the spectrum we want spans the range $4000 \rightarrow 400\text{cm}^{-1}$ (or in near infrared instruments $12,000 \rightarrow 2,500\text{cm}^{-1}$). If the processor is told to generate a cosine function for radiation frequency 4000cm^{-1} (wavelength $2.5\text{ }\mu\text{m}$) it can then search the memorised interferogram for its presence to produce a function F. The search can then be made for radiation of frequency 3999cm^{-1} and

so on down to 400cm^{-1} . The outcome will be the results shown in Figure

A series of points define how much of the various waveforms are present in the interferogram. Put another way, the FT processor acts as a 'frequency analyzer'. If the computer is then told to join up the points the result is a recognisable spectrum.


Figure 8.

[Now go to your machine and run a spectrum of anything you like at 4cm^{-1} resolution. Put it up on the screen and then expand the cm^{-1} axis. Once very well expanded you will see that indeed the spectrum is a series of straight lines joining points].

So – we have our interferometer and hence, as we shall see our interferogram. The next problem is to measure the path difference. A ruler is no good, the path difference has to be measured to an accuracy limited only by the precision of the

engineering in the interferometer itself. The way this is done is to use the interferometer itself to do the job and to feed it with a source of precisely known wavelength. The He/Ne laser has a wavelength near 632.8nm and is known to very high precision. If a He/Ne laser beam passes through the interferometer, its intensity at the output will generate a superb simple cosine with a peak in intensity every 632.8nm of change in the path difference. How this can be done is shown below in Figure 14.


Figure 14. An example of an early FTIR - the 1700 Series Perkin Elmer instrument. J= Jacquinot stop, J₁ is real, J₂ is imaginary. Both mirrors M are fixed

The He/Ne laser generates a parallel beam reflected along the axis of the interferometer by tiny mirror M_{L1}. The beam leaving the interferometer is picked up by mirror M_{L2} and passed to a detector [or rather 2 detectors. If polarised light is used, (thus polariser P) the interfered radiation is split into two polarisation analysed signals, it is possible for the electronics to work out whether the path difference is increasing or decreasing]. The output of the detector(s) is a cosine. As the signal switches from positive to negative about the mean, a trigger circuit is operated providing a pulse twice every 632.8nm of change in path difference. As we shall see, this pulse train is in effect used to measure the path difference.


Figure 15. Output from the He/Ne detector (usually called the fringe detector(s)).

In modern instruments, the He/Ne laser beam is moved out of the optic axis of the interferometer so that the infrared beam is not contaminated with red He/Ne radiation. It is usually passed along an optical bath above or below the infrared beam.

The radiation leaving the interferometer is then passed through the sample area before being focussed onto the detector. Figure 16a and b show the optical diagram and schematic of a complete instrument.


Figure 16a.


Figure 16b. A modern compact FTIR.
X is the pivot axis of the swinging mirrors A & B.

This is all very well, but how is the interferometer actually aligned? The job is done roughly at the factory, but in use it needs very frequent optimisation. This is done by tweaking the orientation of the fixed mirror M_M over tiny angles using electric motors or solenoids driven by software. In fact, in some research grade interferometers several $^{4}\text{He}/^{20}\text{Ne}$ laser beams are used and the alignment is maintained at perfection continuously.

Detectors

Various detectors are encountered in FTIRs, FTNIRs and F-T Raman instruments. The vast majority act as photo resistors i.e. they have a very high resistance in the dark and this falls as light falls on them. The most sensitive are the Ge and InGaAs semi-conductor devices. In the dark, they can have resistance as high as $3 \times 10^8 \Omega$. Measuring resistance at these high values and doing so rapidly requires the use of very clever electronics if spurious signal (noise) is not to be added to the detected signal. All semi-conductor detectors show an 'absorption edge' i.e. they will ignore radiation longer than a characteristic wavelength.

Cooling detectors invariably reduces the amount of noise they develop [we are all familiar with noise. If you listen to your car radio you will be aware that in some areas the sound is clear, in others it is accompanied by hiss – this is "noise" in electronics language. The ratio of signal: noise (S:N) is critical in the

measurement of electrical signals and also in the pleasure in listening to your car radio or mobile phone], but unfortunately cooling shifts the absorption edge towards shorter wavelength. Some detectors give adequate performance (an acceptable useful S:N ratio) at room temperature but others must be cooled.

In an FTIR, one normally finds a TGS or similar detector for ordinary use. If better results [or more often, you need results on samples or sample accessories that transmit very little light] are required, one then resorts to the use of a cryogenically cooled mercury cadmium telluride semi-conductor detector. Thus, these are invariably used in infrared microscopy and very often in diffuse reflection experiments.

The output from the detector goes to a preamplifier where it is converted into a voltage signal varying with time. This signal has to be digitised and the job is usually done with a dedicated 'analogue to digital' converter chip. These devices will measure the input voltage in binary numbers with 8, 16, 20 or 32 or more digits. Clearly, the number of useful digits is governed by the quality (S:N ratio) of the signal. 16 and 20 bit devices are frequently used.

An A/D converter will carry out the conversion only when it is told to do so – it needs a pulse (or 'handshake') to tell it to do the job. The pulse is provided by the He/Ne fringe detector and trigger circuit. So, the A/D converter makes its measurement every 316.4nm of path difference. Clever, ain't it? Thus, the He/Ne makes it possible to measure the interferogram at ultra precise path difference intervals. If the temperature changes, or you lean on the FTIR, the interferometer will very slightly distort and hence the interferogram will shift. However, the He/Ne beam will ALSO shift helping to correct for the distortion.

The digitised interferogram is now fed to the F-T processor. This can be a dedicated chip or a PC. Obviously, a dedicated chip has its advantages (simplicity, lower cost etc) but the PC increases the versatility enormously. The block diagram looks like Figure 16b.

Recording spectra

When you record a spectrum you are, in fact recording two and ratioing them. Let me explain. You empty the sample area (or accessory) and record a BACKGROUND – you then insert the sample, run again and the computer produces a spectrum. What actually happens is that you first record the emission spectrum of source attenuated by the losses in the instrument, the response characteristics of the detector and any absorption along the optical path. The resultant BACKGROUND looks very different from the true emission characteristics of the source. Thus, the beam splitter characteristics and the performance of the detector attenuate the high frequency end of the spectrum. Meanwhile, H₂O+CO₂ in the optical path between the source and the detector cause sharp absorptions [Run the 'BACKGROUND' on your instrument and display it on the screen. You will see intense absorption around 2300cm⁻¹ and 670cm⁻¹ – these are due to CO₂. The forest of bands around 3200 and 1600cm⁻¹ are due to the vibration rotation spectrum of water – for an explanation refer to my piece in [Edition 3, Volume 5](#) earlier this year].

The next spectrum you then record is with the sample in place and again is the emission spectrum of the source, attenuated as before but ALSO by the sample. Usually this is not displayed but rather the second spectrum is immediately ratioed against the first on a point-by-point basis. The result is to show the

difference between the two i.e. the spectrum of the sample which is what you want.

Once you understand what happens several sources of possible error become clear –


You can't get meaningful results if you ratio nothing against nothing! So, if the background emission is non-existent, as it is over the CO₂ absorption near 2300cm⁻¹, you get no useful results.

If you record the background under one set of conditions (resolution, number of scans etc) and the spectrum under another, the ratioing will be meaningless. It turns out that if the number of scans in the BACKGROUND is greater than that in the sample spectrum you are OK. You can, however, change the conditions between BACKGROUND and sample spectrum unwittingly.


If you breath into the sample area as you insert the sample you will increase the CO₂ and water vapour in the optical path. As a result, the BACKGROUND you are using is not appropriate and spurious bands will appear in your spectrum. If your instrument is purged with dry nitrogen, opening the sample area will introduce CO₂ + H₂O so the BACKGROUND and sample spectrum must again contain unwanted differences. This article is getting a little long so let me come back to these matters in another piece in a future edition.

Resolution

Like S:N ratio resolution as actually visualised by the user is somewhat subjective. We all have a sort of feel for RESOLUTION – it's the ability to pick out one band from another lying close to it. So, we want to be sure that if two bands lie xcm⁻¹ apart we can see them i.e.

Lord Rayleigh way back in the 1800's defined resolution in terms of the appearance of the spectra but now that we have incredibly monochromatic laser sources people define resolution in terms of 'bandwidth'. If you take a truly monochromatic source and pass its emission through an instrument you will actually see

The width of the band at exactly $\frac{1}{2}$ of its height is defined as the resolution. The same criteria apply to absorptions – if the absorption is very sharp and almost monochromatic e.g. a vibration rotation band of a gas at a low pressure, the $W_{1/2}$ value of the absorption band will be its measured resolution.

The instrument makers offer us resolution values on FTIRs. Do the offered values match what I have said above? Well, the answer is yes and no! The industry accepts a rather different set of criteria to define 'resolution'. The standard is that the 'resolution' is related to the data point density in the Fourier Transformed Spectrum, the standard being that the resolution is $\frac{1}{2}$ the data point density. Let me give you an example --if you set the resolution at 2cm^{-1} , the instrument will compute a spectral data point every wavenumber 4000, 3999, 3998.....401, 400cm^{-1} .

The data points will then be joined up and that's what you get. The quality of spectrum you actually see depends on the wavenumbers where the data points are calculated and the wavenumber of the emission or absorption bands but the industry standard is reasonable.

One obvious question is – why not record all spectra at the best resolution of which my equipment is capable? The answer is time (or more strictly S:N ratio).


Remember, the resolution is standardised as a half of the data point density – so the number of data points between 4400 and 400cm^{-1} will be a modest 2000 at 4cm^{-1} resolution, but will rise to 16,000 at 0.5cm^{-1} resolution. To be able to calculate more data points, more points must be measured i.e. the interferogram must be longer. In fact, the way it works out, the length of the interferogram – the distance scanned by the moving mirror equals $\frac{1}{\text{resolution}}$. So – a machine capable of 0.1cm^{-1} resolution requires that M_M will move by at least 10cm, a long way remembering the precision of movement required. This is why small cheap machine tend to offer only restricted resolution.

Recording more data points is all very well, but more time is then needed to calculate the points in the spectrum. To sum up – using a better resolution than is required is simply not worth it because every measurement costs time.

Apodization Functions

When your FTIR scans, it actually records an interferogram, which is subsequently digitised, and then Fourier Transformed into a spectrum ([See Article 1 for details](#)).

The form of the interferogram is

The vertical axis is never labelled i.e. it is in arbitrary units BUT in reality it is the output voltage or current of the detector. Unlike a spectrum, which shows peaks against a background, the interferogram shows swings of output about a non-zero mean. At the end of each scan, the machine reverses its direction. Inevitably the reversal means that the scan must be slowed down, stopped, reversed and then sped up again. As the optical delay is slowed, stopped and then accelerated, the output data is corrupted, so the instrument is told to switch off the detector until the reversal is complete and the movement has become nice and smooth. So to put it diagrammatically we have -


Figure 2.

And when you display the interferogram you only 'see' the good bit, the span between the two switching intervals.

Fine – but there is a problem. If you switch off a signal, it will naturally fall rapidly to zero i.e. there will be a 'step' in the output value of the detector circuit. Unfortunately, the Fourier Transform processor will process this step and will make spectral nonsense out of it. The trick is to minimise this nuisance by lowering the output along a chosen curve – an "apodization function". Each function has a different shape eg.,


Figure 3.

To demonstrate the effect of these, we recorded spectra of polystyrene under various conditions and compared them. All spectra were recorded with a TGS detector at 4cm^{-1} using 4 co-added scans. The backgrounds and the spectra were each recorded under identical conditions.

	Spectrum 1	Strong (Beer Norton) Apodization	See Figure 1
Compare	Spectrum 2	Weak apodization	See Figure 4a and 4b
	Spectrum 3	Boxcar – no apodization	See Figure 5a and 5b


Figure 4. (a) background (b) PS film both with a weak

apodization


Figure 5. (a) background (b) PS film both with a boxcar apodization

The spectra of polystyrene showed little or no effect as one changes the apodization however, where bands are sharp one can see effects. Look at the backgrounds and you will see that reducing or switching off the apodization can cause the sharp vibration rotation bands of water vapour to alter in apparent position and shape. So – if you are going to record spectra of gases especially at higher resolutions or of crystalline solids with very sharp bands you must be careful about apodization. If too low the bands will appear to be better resolved but 'ringing' can occur.


Figure 6.

In most analytical spectroscopy there is little problem. Since the bands are broad, the effect of using high apodization is not disastrous and it is safer to use it rather than not.

Smoothing

Older spectrometers offered simple smoothing procedures eg., the instrument took 3, 5 or even more spectral data points, averaged them and plotted out the average. It then moved by one cm^{-1} increment (usually $\frac{1}{2}$ resolution value) and repeated. These days more sophisticated functions are used but the effect is similar – you reduce the noise by averaging but unfortunately you lose resolution.

Let's try an experiment. Take a piece of Al foil (kitchen quality is just fine) and punch a hole in it 1 to 1.5mm in diameter with a needle or sharp pencil. Place this over your polystyrene film and put the sandwich in the spectrometer after you have run a background. Make sure the hole is roughly in the middle of the IR beam. Run a spectrum. You should see something like this –


Figure 7.

That is the whole spectrum is squashed up – very little light passes the hole in the foil. Now expand the spectrum to fill the screen – Figure 8.


Figure 8.

The spectrum is fine but the result is noisy. This problem can arise where for some reason the energy is very restricted (examples include the use of some spectral reflection techniques, spectra on samples held in some high or low temperature cells or in IR microscopy). The correct way to deal with it is to co-add many scans [Remember, the S:N ratio – improves as $\sqrt{\text{v\# of scans}}$]. Either this may not be possible or just take too long, or the experiment is finished and you want to clean up the results you have and not do the job again. The application of smoothing is worth considering.

In Figure 9 we show the effect of smoothing the data in Figure 8.


Figure 9.

There is a loss of resolution but the noise is indeed reduced. In this case the loss of resolution is acceptable but in some other cases it may not be.

Absorbance vs Transmittance

In our experience most analytical users of infrared record spectra in transmission (0-100% transmission) and never use absorbance. The theory behind this is –

$$Tr = \frac{I_{Tr}}{I_0}, \quad \%Tr = 100 \cdot \frac{I_{Tr}}{I_0}$$

$$Ab = -\log_{10}\left(\frac{I_{Tr}}{I_0}\right)$$

The negative sign means that the peaks in Absorbance are the negative ones in transmission.

In Figure 10 we show part of the spectrum of polystyrene presented in transmission (bands downwards) and absorbance (bands upwards).



Figure 10.

Several points are clear –

- Absorbance spectra look much better than transmittance. The bands appear to be sharper and clearer.
- The width at half height $W_{1/2}$ is narrow in transmission.

The band at 1451cm^{-1} in polystyrene has $W_{1/2}$ of $\sim 18\text{cm}^{-1}$ in transmittance and $\sim 2\frac{1}{2}\text{cm}^{-1}$ in absorbance. If you check the vibration rotation gas phase spectrum of water vapour in the background you will find less of an effect eg., the band at 1050cm^{-1} narrowed from $W_{1/2} = 6.5$ to about 5cm^{-1} .

So the question must arise – why not use absorbance more often as a method of presentation.

You will note that in Figure 10b, the absorbance scale goes from 0 to 7 i.e. 100 to about 0.1%.

REF: P.J. Hendra & F. Birembaut, *Int.J.Vibr.Spec.*, [www.ijvs.com] **6**, 2, 3 (2002)

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